



SAIKANG Medical

Healthier your life

STOCK CODE: 870098

SKE013-10

LAB CHAIR



SKE013-10

FEATURES AND HIGHLIGHTS

Technical parameters

▪ Stool surface diameter	350mm
▪ Thickness	50mm
▪ Including wheel height	565-805mm
▪ Back high	750mm



Structural design

Chrome plated structure, corrosion resistance, no rust, no water, no moisture, safe and reliable.



Stool surface

PU leather, polyurethane foam seat, soft and comfortable.



Gas Spring

Stainless steel material, $\Phi 450$ mm, it allows the feet to be placed more comfortably and sitting more comfortably.



Chassis

Stainless steel material, $\Phi 450$ mm, it allows the feet to be placed more comfortably and sitting more comfortably.



Castors

Five nylon castors, convenient, wear-resistant, sliding silent.



Function

The seat cushion can rotate around and lift up and down freely. The chassis can adjust the foot height up and down.





CONSTANT INNOVATION



Design



Research



Manufacturing



Maintenance



Service



WORLDWIDE BRANCHES

South America Branch

Saikang Medical Technology, Inc

W3400 Corrientes, Argentina
TEL: +55 (11) 2758-0438
Rua Universal, Hollywood,
São Bernardo do Campo,
Brazil
TEL: +55 (11) 2758-0438

East Africa Branch

Saikang Medical Technology Ltd

Mezzanine Floor, Taj Tower,
upper hill, nairobi, Kenya
TEL: +254714084806

South East Asia Branch

Saikang Medical Technology Ltd

Sala Thammasop Subdistrict
Thawi
Watthana District
Bangkok Thailand
TEL: +66 642741046

US Branch

Healthward Technology, INC

17800 CASTLETON ST STE
180,
CITY OF INDUSTRY CA 91748
TEL: +1 (603)416-3063

Russia Branch

LLC 'SAIKANG'

454004, Russian Federation,
Chelyabinsk,
Universitetskaya
Naberezhnaya, Russia
TEL: +7 (919)300-00-06

Mid East Branch

Saikang Medical Supplies (Mid East) LLC

Musaffah
P.O Box 53200, Abu Dhabi,
UAE
Tel: +971 585710652

SAIKANG Medical

Jiangsu Saikang Medical Equipment Co., Ltd.
STOCK CODE: 870098

📍 No. 35, Lehong Road · Zhangjiagang (Jiangsu) · China

☎ +86 18021231901 · 🕒 24/7

✉ export@saikangmedical.com

🌐 www.saikangmedical.com



Company video



General catalogue